


Genotyping Report

Strain ID	T024709	Strain Type	CKO(Cas9)	Genetic Background	C57BL/6J
Designer	Ya'nan Xu	Gene Name	<i>Med24</i>		

1. Strategy of Genotyping




Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.
 Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.
 Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.
 Note: The sizes of WT and Targeted band are shown below.

2. Primer Information

PCR No.	Primer Name	Primer Name	Sequence	Band Size
①(5'arm)	F1	JS25003-Med24-5wt-tF1	GCCCCGGCACAGTTAATCTTGATAA	WT: 258bp Targeted:362bp
	R1	JS25003-Med24-5wt-tR1	AAAGAATCACCCAGAATGACGGTC	
②(3'arm)	F2	JS25003-Med24-3wt-tF1	TAAGGGTTCTCCATTGACAGGAATG	WT: 306bp Targeted: 407bp
	R2	JS25003-Med24-3wt-tR1	GCCAATTAACCTCCATGTAGCACAG	

3. Gel Image & Conclusion



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component	Volume (μl)	
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)	12.5	
2	ddH ₂ O	9.5	
3	Primer A(10pmol/μl)	1	
4	Primer B(10pmol/μl)	1	
5	Template(≈100ng/μl)	1	
PCR program			
Seg.	Temp.	Time	Cycle
1	95°C	5min	
2	98°C	30s	20×
3	65°C (-0.5°C/cycle)	30s	
4	72°C	45s	
5	98°C	30s	20×
6	55°C	30s	
7	72°C	45s	
8	72°C	5min	
9	10°C	hold	

Reviewer:Ting Sun

Date:2020-4-7