


Genotyping Report

Strain ID	T005122	Strain Type	CKO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Ya'nan Xu	Gene Name	<i>Usp51</i>		

1. Strategy of Genotyping



2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
①(5'arm)	T005122-F1	GAGTAGAGGCTGGAAGCATAGGTAA	WT:251bp Targeted:308bp
	T005122-R1	ACTGTCCTGAGCTTTAGGGTGAA	
②(3'arm)	T005122-F2	CTTATGTCAATGCAGATAACTATATCG	WT: 324bp Targeted:381bp
	T005122-R2	GATAAGGAATTAGCAGGACCGTTA	

3. Gel Image & Conclusion



4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus , (Vazyme P112-03)		12.5
2	ddH ₂ O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95 °C	5min	
2	98 °C	30s	20×
3	65 °C (-0.5 °C/cycle)	30s	
4	72 °C	45s	
5	98 °C	30s	20×
6	55 °C	30s	
7	72 °C	45s	
8	72 °C	5min	
9	10 °C	hold	

Reviewer:Ting Sun

Date:2021-6-9