


## Genotyping Report

Strain ID	T019281	Strain Type	CKO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Xiaomin Shi	Gene Name	<i>Cbln1</i>		

### 1. Strategy of Genotyping



Wild type: ①PCR reaction obtains a single WT band; ②PCR reaction obtains a single WT band.  
 Heterozygote: ①PCR reaction obtains a WT band and a Targeted band; ②PCR reaction obtains a WT band and a Targeted band.  
 Homozygote: ①PCR reaction obtains a single Targeted band; ②PCR reaction obtains a single Targeted band.  
 Note: The sizes of WT and Targeted band are shown below.

### 2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
① (5'arm)	T019281-F1	GAATACTGCCTTTCCTGGACCTG	WT: 353bp Targeted: 458bp
	T019281-R1	GCTTAACTTGTTTCGGCCTTCTG	
② (3'arm)	T019281-F2	GCCTCATAATTGCTCTTTGGC	WT: 239bp Targeted: 345bp
	T019281-R2	GAGCCTGATTGTTGAGCCTAGAAAG	

### 3. Gel Image & Conclusion

--



#### 4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus, (Vazyme P112-03)		12.5
2	ddH <sub>2</sub> O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95°C	5min	
2	98°C	30s	20×
3	65°C (-0.5°C/cycle)	30s	
4	72°C	45s	
5	98°C	30s	20×
6	55°C	30s	
7	72°C	45s	
8	72°C	5min	
9	10°C	hold	