


## Genotyping Report

Strain ID	T030249	Strain Type	KO(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Zifan Lin	Gene Name	<i>Sbf2</i>		


### 1. Strategy of Genotyping



### 2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
PCR①	T030249-F1	GTCTGGAACATTCCATTCTTGGT	WT: 6505bp KO:404bp
	T030249-R1	CTTTCACAGTAGGGCAGGACAGTG	
PCR②	T030249-F2	GATCATCAGGTTACACGGTGTG	WT: 397bp KO:0bp
	T030249-R2	CCAGCCAGAGGGATTAAAGGTGTG	

### 3. Gel Image



#### 4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus , (Vazyme P112-03)		12.5
2	ddH <sub>2</sub> O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95℃	5min	
2	98℃	30s	20×
3	65℃ (-0.5℃/cycle)	30s	
4	72℃	45s	
5	98℃	30s	20×
6	55℃	30s	
7	72℃	45s	
8	72℃	5min	
9	10℃	hold	

Reviewer:Ting Sun

Date:2021-6-16