


Genotyping Report

Strain ID	T050765	Strain Type	KI(Cas9)	Genetic Background	C57BL/6JGpt
Designer	Dongdong Zhang	Gene Name	<i>Trpv1-IRES-EGFP</i>		


1. Strategy of Genotyping



2. Primer Information

PCR No.	Primer No.	Sequence	Band Size
①5'arm	T050765-F1A	ACCATGCTGCCTTTCAGTTTCAG	WT:0bp Targeted:450bp
	T050765-R1A	TTCCAGAGGAACTGCTTCCTTCAC	
②WT	T050765-F2	ACCATGCTGCCTTTCAGTTTCAG	WT:608bp Targeted:1958bp
	T050765-R2	TCCCAGGTAAGATAGTCAGCAGAGC	

3. Gel Image & Conclusion



Note: P:Positive control; WT: Wildtype control; B: Blank control (ddH₂O); M: DNA Ladder

4. PCR Condition

PCR Reaction Component			
Seg.	reaction component		Volume (μl)
1	2 × Taq Master Mix , Dye Plus' (Vazyme P112-03)		12.5
2	ddH ₂ O		9.5
3	Primer A(10pmol/μl)		1
4	Primer B(10pmol/μl)		1
5	Template(≈100ng/μl)		1
PCR program			
Seg.	Temp.	Time	Cycle
1	95°C	5min	20×
2	98°C	30s	
3	65°C (-0.5°C/cycle)	30s	
4	72°C	45s	
5	98°C	30s	20×
6	55°C	30s	
7	72°C	45s	
8	72°C	5min	
9	10°C	hold	